

RUSH CHOPART

Kategorivelger i forhold til aktivitetsnivå, brukervekt og størrelse																				
Vekt	0 - 32 kg		33 - 48 kg		49 - 64 kg		64 - 79 kg		80 - 95 kg		96 - 111 kg		112 - 127 kg		128 - 143 kg		144 - 159 kg		160 - 166 kg	
	Aktivitetsnivå	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact	
Størrelse	0	1	1	2	2	3	3	4	4	5	5	-	-	-	-	-	-	-	-	
22 cm	0	1	1	2	2	3	3	4	4	5	5	-	-	-	-	-	-	-	-	
23 cm	0	1	1	2	2	3	3	4	4	5	5	-	-	-	-	-	-	-	-	
24 cm	-	-	1	2	2	3	3	4	4	5	5	6	6	-	-	-	-	-	-	
25 cm	-	-	-	-	2	3	3	4	4	5	5	6	6	-	-	-	-	-	-	
26 cm	-	-	-	-	2	3	3	4	4	5	5	6	6	7	7	-	-	-	-	
27 cm	-	-	-	-	-	-	3	4	4	5	5	6	6	7	7	-	-	-	-	
28 cm	-	-	-	-	-	-	3	4	4	5	5	6	6	7	7	8	8	-	-	
29 cm	-	-	-	-	-	-	3	4	4	5	5	6	6	7	7	8	8	-	-	
30 cm	-	-	-	-	-	-	-	-	4	5	5	6	6	7	7	8	8	9	9	

Kategorivelger i forhold til aktivitetsnivå, brukervekt og størrelse																	
Vekt	0 - 14 kg		15 - 20 kg		21 - 27 kg		28 - 34 kg		35 - 41 kg		42 - 48 kg		49 - 54 kg		55 - 61 kg		
	Aktivitetsnivå	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact
Størrelse	1P	2P	2P	3P	3P	4P	4P	5P	5P	6P	6P	-	-	-	-	-	
19 cm	1P	2P	2P	3P	3P	4P	4P	5P	5P	6P	6P	-	-	-	-	-	
20 cm	-	-	2P	3P	3P	4P	4P	5P	5P	6P	6P	7P	7P	-	-	-	
21 cm	-	-	2P	3P	3P	4P	4P	5P	5P	6P	6P	7P	7P	8P	8P	-	

Art.nr.	Modell	Lev.art.nr.
759701910	SH Chopart 19cm Cat 1	CHO-19-1-RU
759701920	SH Chopart 19cm Cat 2	CHO-19-2-RU
759701930	SH Chopart 19cm Cat 3	CHO-19-3-RU
759701940	SH Chopart 19cm Cat 4	CHO-19-4-RU
759701950	SH Chopart 19cm Cat 5	CHO-19-5-RU
759701960	SH Chopart 19cm Cat 6	CHO-19-6-RU
759702020	SH Chopart 20cm Cat 2	CHO-20-2-RU
759702030	SH Chopart 20cm Cat 3	CHO-20-3-RU
759702040	SH Chopart 20cm Cat 4	CHO-20-4-RU
759702050	SH Chopart 20cm Cat 5	CHO-20-5-RU
759702060	SH Chopart 20cm Cat 6	CHO-20-6-RU
759702070	SH Chopart 20cm Cat 7	CHO-20-7-RU
759702120	SH Chopart 21cm Cat 2	CHO-21-2-RU
759702130	SH Chopart 21cm Cat 3	CHO-21-3-RU
759702140	SH Chopart 21cm Cat 4	CHO-21-4-RU
759702150	SH Chopart 21cm Cat 5	CHO-21-5-RU
759702160	SH Chopart 21cm Cat 6	CHO-21-6-RU
759702170	SH Chopart 21cm Cat 7	CHO-21-7-RU
759702180	SH Chopart 21cm Cat 8	CHO-21-8-RU
759702200	SH Chopart 22cm Cat 0	CHO-22-0-RU
759702210	SH Chopart 22cm Cat 1	CHO-22-1-RU
759702220	SH Chopart 22cm Cat 2	CHO-22-2-RU
759702230	SH Chopart 22cm Cat 3	CHO-22-3-RU
759702240	SH Chopart 22cm Cat 4	CHO-22-4-RU
759702250	SH Chopart 22cm Cat 5	CHO-22-5-RU
759702300	SH Chopart 23cm Cat 0	CHO-23-0-RU
759702310	SH Chopart 23cm Cat 1	CHO-23-1-RU
759702320	SH Chopart 23cm Cat 2	CHO-23-2-RU
759702330	SH Chopart 23cm Cat 3	CHO-23-3-RU

759702340	SH Chopart 23cm Cat 4	CHO-23-4-RU
759702350	SH Chopart 23cm Cat 5	CHO-23-5-RU
759702410	SH Chopart 24cm Cat 1	CHO-24-1-RU
759702420	SH Chopart 24cm Cat 2	CHO-24-2-RU
759702430	SH Chopart 24cm Cat 3	CHO-24-3-RU
759702440	SH Chopart 24cm Cat 4	CHO-24-4-RU
759702450	SH Chopart 24cm Cat 5	CHO-24-5-RU
759702460	SH Chopart 24cm Cat 6	CHO-24-6-RU
759702520	SH Chopart 25cm Cat 2	CHO-25-2-RU
759702530	SH Chopart 25cm Cat 3	CHO-25-3-RU
759702540	SH Chopart 25cm Cat 4	CHO-25-4-RU
759702550	SH Chopart 25cm Cat 5	CHO-25-5-RU
759702560	SH Chopart 25cm Cat 6	CHO-25-6-RU
759702620	SH Chopart 26cm Cat 2	CHO-26-2-RU
759702630	SH Chopart 26cm Cat 3	CHO-26-3-RU
759702640	SH Chopart 26cm Cat 4	CHO-26-4-RU
759702650	SH Chopart 26cm Cat 5	CHO-26-5-RU
759702660	SH Chopart 26cm Cat 6	CHO-26-6-RU
759702670	SH Chopart 26cm Cat 7	CHO-26-7-RU
759702730	SH Chopart 27cm Cat 3	CHO-27-3-RU
759702740	SH Chopart 27cm Cat 4	CHO-27-4-RU
759702750	SH Chopart 27cm Cat 5	CHO-27-5-RU
759702760	SH Chopart 27cm Cat 6	CHO-27-6-RU
759702770	SH Chopart 27cm Cat 7	CHO-27-7-RU
759702830	SH Chopart 28cm Cat 3	CHO-28-3-RU
759702840	SH Chopart 28cm Cat 4	CHO-28-4-RU
759702850	SH Chopart 28cm Cat 5	CHO-28-5-RU
759702860	SH Chopart 28cm Cat 6	CHO-28-6-RU
759702870	SH Chopart 28cm Cat 7	CHO-28-7-RU
759702880	SH Chopart 28cm Cat 8	CHO-28-8-RU
759702930	SH Chopart 29cm Cat 3	CHO-29-3-RU
759702940	SH Chopart 29cm Cat 4	CHO-29-4-RU
759702950	SH Chopart 29cm Cat 5	CHO-29-5-RU
759702960	SH Chopart 29cm Cat 6	CHO-29-6-RU
759702970	SH Chopart 29cm Cat 7	CHO-29-7-RU
759702980	SH Chopart 29cm Cat 8	CHO-29-8-RU
759703040	SH Chopart 30cm Cat 4	CHO-30-4-RU
759703050	SH Chopart 30cm Cat 5	CHO-30-5-RU
759703060	SH Chopart 30cm Cat 6	CHO-30-6-RU
759703070	SH Chopart 30cm Cat 7	CHO-30-7-RU
759703080	SH Chopart 30cm Cat 8	CHO-30-8-RU
759703090	SH Chopart 30cm Cat 9	CHO-30-9-RU