

RUSH RAMPAGE

Kategorivelger i forhold til aktivitetsnivå, brukervekt og størrelse																		
Vekt	0 - 48 kg		49 - 64 kg		64 - 79 kg		80 - 95 kg		96 - 111 kg		112 - 127 kg		128 - 143 kg		144 - 159 kg		160 - 166 kg	
	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact
22 cm	1	2	2	3	3	4	4	5	5	6	6	-	-	-	-	-	-	-
23 cm	1	2	2	3	3	4	4	5	5	6	6	-	-	-	-	-	-	-
24 cm	-	-	2	3	3	4	4	5	5	6	6	-	-	-	-	-	-	-
25 cm	-	-	2	3	3	4	4	5	5	6	6	7	7	-	-	-	-	-
26 cm	-	-	2	3	3	4	4	5	5	6	6	7	7	8	8	-	-	-
27 cm	-	-	-	-	3	4	4	5	5	6	6	7	7	8	8	9	9	-
28 cm	-	-	-	-	3	4	4	5	5	6	6	7	7	8	8	9	9	-
29 cm	-	-	-	-	3	4	4	5	5	6	6	7	7	8	8	9	9	-

Art.nr.	Modell	Lev.art.nr.
759122100	RUSH RAMPAGE 22cm cat 1	RAM-22-1-RU
759122200	RUSH RAMPAGE 22cm cat 2	RAM-22-2-RU
759122300	RUSH RAMPAGE 22cm cat 3	RAM-22-3-RU
759122400	RUSH RAMPAGE 22cm cat 4	RAM-22-4-RU
759122500	RUSH RAMPAGE 22cm cat 5	RAM-22-5-RU
759123100	RUSH RAMPAGE 23cm cat 1	RAM-23-1-RU
759123200	RUSH RAMPAGE 23cm cat 2	RAM-23-2-RU
759123300	RUSH RAMPAGE 23cm cat 3	RAM-23-3-RU
759123400	RUSH RAMPAGE 23cm cat 4	RAM-23-4-RU
759123500	RUSH RAMPAGE 23cm cat 5	RAM-23-5-RU
759123600	RUSH RAMPAGE 23cm cat 6	RAM-23-6-RU
759124100	RUSH RAMPAGE 24cm cat 1	RAM-24-1-RU
759124200	RUSH RAMPAGE 24cm cat 2	RAM-24-2-RU
759124300	RUSH RAMPAGE 24cm cat 3	RAM-24-3-RU
759124400	RUSH RAMPAGE 24cm cat 4	RAM-24-4-RU
759124500	RUSH RAMPAGE 24cm cat 5	RAM-24-5-RU
759124600	RUSH RAMPAGE 24cm cat 6	RAM-24-6-RU
759125200	RUSH RAMPAGE 25cm cat 2	RAM-25-2-RU
759125300	RUSH RAMPAGE 25cm cat 3	RAM-25-3-RU
759125400	RUSH RAMPAGE 25cm cat 4	RAM-25-4-RU
759125500	RUSH RAMPAGE 25cm cat 5	RAM-25-5-RU
759125600	RUSH RAMPAGE 25cm cat 6	RAM-25-6-RU
759125700	RUSH RAMPAGE 25cm cat 7	RAM-25-7-RU
759126200	RUSH RAMPAGE 26cm cat 2	RAM-26-2-RU
759126300	RUSH RAMPAGE 26cm cat 3	RAM-26-3-RU
759126400	RUSH RAMPAGE 26cm cat 4	RAM-26-4-RU
759126500	RUSH RAMPAGE 26cm cat 5	RAM-26-5-RU
759126600	RUSH RAMPAGE 26cm cat 6	RAM-26-6-RU
759126700	RUSH RAMPAGE 26cm cat 7	RAM-26-7-RU
759127300	RUSH RAMPAGE 27cm cat 3	RAM-27-3-RU
759127400	RUSH RAMPAGE 27cm cat 4	RAM-27-4-RU
759127500	RUSH RAMPAGE 27cm cat 5	RAM-27-5-RU
759127600	RUSH RAMPAGE 27cm cat 6	RAM-27-6-RU
759127700	RUSH RAMPAGE 27cm cat 7	RAM-27-7-RU
759128300	RUSH RAMPAGE 28cm cat 3	RAM-28-3-RU
759128400	RUSH RAMPAGE 28cm cat 4	RAM-28-4-RU
759128500	RUSH RAMPAGE 28cm cat 5	RAM-28-5-RU
759128600	RUSH RAMPAGE 28cm cat 6	RAM-28-6-RU
759128700	RUSH RAMPAGE 28cm cat 7	RAM-28-7-RU
759128800	RUSH RAMPAGE 28cm cat 8	RAM-28-8-RU
759128900	RUSH RAMPAGE 28cm cat 9	RAM-28-9-RU
759129300	RUSH RAMPAGE 29cm cat 3	RAM-29-3-RU
759129400	RUSH RAMPAGE 29cm cat 4	RAM-29-4-RU
759129500	RUSH RAMPAGE 29cm cat 5	RAM-29-5-RU
759129600	RUSH RAMPAGE 29cm cat 6	RAM-29-6-RU
759129700	RUSH RAMPAGE 29cm cat 7	RAM-29-7-RU
759129800	RUSH RAMPAGE 29cm cat 8	RAM-29-8-RU
759129900	RUSH RAMPAGE 29cm cat 9	RAM-29-9-RU