

RUSH ROVER

Kategoriøvelser i forhold til aktivitetsnivå, brukervekt og størrelse																		
Vekt	0 - 48 kg		49 - 64 kg		64 - 79 kg		80 - 95 kg		96 - 111 kg		112 - 127 kg		128 - 143 kg		144 - 159 kg		160 - 166 kg	
	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact
22 cm	1	2	2	3	3	4	4	4	-	-	-	-	-	-	-	-	-	-
23 cm	1	2	2	3	3	4	4	4	5	5	-	-	-	-	-	-	-	-
24 cm	1	2	2	3	3	4	4	4	5	5	-	-	-	-	-	-	-	-
25 cm	1	2	2	3	3	4	4	4	5	5	6	6	-	-	-	-	-	-
26 cm	-	-	2	3	3	4	4	4	5	5	6	6	7	7	-	-	-	-
27 cm	-	-	2	3	3	4	4	4	5	5	6	6	7	7	-	-	-	-
28 cm	-	-	-	-	3	4	4	4	5	5	6	6	7	7	8	8	-	-
29 cm	-	-	-	-	3	4	4	4	5	5	6	6	7	7	8	8	-	-
30 cm	-	-	-	-	-	-	-	4	5	5	6	6	7	7	8	8	9	9

Art.nr.	Modell	Lev.art.nr.
759522100	RUSH ROVER 22cm Cat 1	ROV-22-1-RU
759522200	RUSH ROVER 22cm Cat 2	ROV-22-2-RU
759522300	RUSH ROVER 22cm Cat 3	ROV-22-3-RU
759522400	RUSH ROVER 22cm Cat 4	ROV-22-4-RU
759523100	RUSH ROVER 23cm Cat 1	ROV-23-1-RU
759523200	RUSH ROVER 23cm Cat 2	ROV-23-2-RU
759523300	RUSH ROVER 23cm Cat 3	ROV-23-3-RU
759523400	RUSH ROVER 23cm Cat 4	ROV-23-4-RU
759523500	RUSH ROVER 23cm Cat 5	ROV-23-5-RU
759524100	RUSH ROVER 24cm Cat 1	ROV-24-1-RU
759524200	RUSH ROVER 24cm Cat 2	ROV-24-2-RU
759524300	RUSH ROVER 24cm Cat 3	ROV-24-3-RU
759524400	RUSH ROVER 24cm Cat 4	ROV-24-4-RU
759524500	RUSH ROVER 24cm Cat 5	ROV-24-5-RU
759525100	RUSH ROVER 25cm Cat 1	ROV-25-1-RU
759525200	RUSH ROVER 25cm Cat 2	ROV-25-2-RU
759525300	RUSH ROVER 25cm Cat 3	ROV-25-3-RU
759525400	RUSH ROVER 25cm Cat 4	ROV-25-4-RU
759525500	RUSH ROVER 25cm Cat 5	ROV-25-5-RU
759525600	RUSH ROVER 25cm Cat 6	ROV-25-6-RU
759526200	RUSH ROVER 26cm Cat 2	ROV-26-2-RU
759526300	RUSH ROVER 26cm Cat 3	ROV-26-3-RU
759526400	RUSH ROVER 26cm Cat 4	ROV-26-4-RU
759526500	RUSH ROVER 26cm Cat 5	ROV-26-5-RU
759526600	RUSH ROVER 26cm Cat 6	ROV-26-6-RU
759526700	RUSH ROVER 26cm Cat 7	ROV-26-7-RU
759527200	RUSH ROVER 27cm Cat 2	ROV-27-2-RU
759527300	RUSH ROVER 27cm Cat 3	ROV-27-3-RU
759527400	RUSH ROVER 27cm Cat 4	ROV-27-4-RU
759527500	RUSH ROVER 27cm Cat 5	ROV-27-5-RU
759527600	RUSH ROVER 27cm Cat 6	ROV-27-6-RU
759527700	RUSH ROVER 27cm Cat 7	ROV-27-7-RU
759528300	RUSH ROVER 28cm Cat 3	ROV-28-3-RU
759528400	RUSH ROVER 28cm Cat 4	ROV-28-4-RU
759528500	RUSH ROVER 28cm Cat 5	ROV-28-5-RU
759528600	RUSH ROVER 28cm Cat 6	ROV-28-6-RU
759528700	RUSH ROVER 28cm Cat 7	ROV-28-7-RU
759528800	RUSH ROVER 28cm Cat 8	ROV-28-8-RU
759529300	RUSH ROVER 29cm Cat 3	ROV-29-3-RU
759529400	RUSH ROVER 29cm Cat 4	ROV-29-4-RU
759529500	RUSH ROVER 29cm Cat 5	ROV-29-5-RU
759529600	RUSH ROVER 29cm Cat 6	ROV-29-6-RU
759529700	RUSH ROVER 29cm Cat 7	ROV-29-7-RU
759529800	RUSH ROVER 29cm Cat 8	ROV-29-8-RU
759530400	RUSH ROVER 30cm Cat 4	ROV-30-4-RU
759530500	RUSH ROVER 30cm Cat 5	ROV-30-5-RU
759530600	RUSH ROVER 30cm Cat 6	ROV-30-6-RU
759530700	RUSH ROVER 30cm Cat 7	ROV-30-7-RU
759530800	RUSH ROVER 30cm Cat 8	ROV-30-8-RU
759530900	RUSH ROVER 30cm Cat 9	ROV-30-9-RU