

# EASY RUN

Kategorivelger i forhold til brukervekt og aktivitetsnivå							
Aktivitet / Vekt	44 - 49 kg	50 - 59 kg	60 - 69 kg	70 - 79 kg	80 - 89 kg	90 - 99 kg	100 - 110 kg
Moderate impact		1	2	3	4	5	6
High impact	1	2	3	4	5	6	-

Art.nr.	Modell	Lev.art.nr.
722500110	EASY RUN CAT. 1 - WATER	1A920-1E
722500120	EASY RUN CAT. 1 - TAG	1A920-1G
722500130	EASY RUN CAT. 1 - BLACK	1A920-1N
722500140	EASY RUN CAT. 1 - ORANGE	1A920-1O
722500150	EASY RUN CAT. 1 - USA	1A920-1U
722500160	EASY RUN CAT. 1 - PURPLE	1A920-1V
722500210	EASY RUN CAT. 2 - WATER	1A920-2E
722500220	EASY RUN CAT. 2 - TAG	1A920-2G
722500230	EASY RUN CAT. 2 - BLACK	1A920-2N
722500240	EASY RUN CAT. 2 - ORANGE	1A920-2O
722500250	EASY RUN CAT. 2 - USA	1A920-2U
722500260	EASY RUN CAT. 2 - PURPLE	1A920-2V
722500310	EASY RUN CAT. 3 - WATER	1A920-3E
722500320	EASY RUN CAT. 3 - TAG	1A920-3G
722500330	EASY RUN CAT. 3 - BLACK	1A920-3N
722500340	EASY RUN CAT. 3 - ORANGE	1A920-3O
722500350	EASY RUN CAT. 3 - USA	1A920-3U
722500360	EASY RUN CAT. 3 - PURPLE	1A920-3V
722500410	EASY RUN CAT. 4 - WATER	1A920-4E
722500420	EASY RUN CAT. 4 - TAG	1A920-4G
722500430	EASY RUN CAT. 4 - BLACK	1A920-4N
722500440	EASY RUN CAT. 4 - ORANGE	1A920-4O
722500450	EASY RUN CAT. 4 - USA	1A920-4U
722500460	EASY RUN CAT. 4 - PURPLE	1A920-4V
722500510	EASY RUN CAT. 5 - WATER	1A920-5E
722500520	EASY RUN CAT. 5 - TAG	1A920-5G
722500530	EASY RUN CAT. 5 - BLACK	1A920-5N
722500540	EASY RUN CAT. 5 - ORANGE	1A920-5O
722500550	EASY RUN CAT. 5 - USA	1A920-5U
722500560	EASY RUN CAT. 5 - PURPLE	1A920-5V
722500610	EASY RUN CAT. 6 - WATER	1A920-6E
722500620	EASY RUN CAT. 6 - TAG	1A920-6G
722500630	EASY RUN CAT. 6 - BLACK	1A920-6N
722500640	EASY RUN CAT. 6 - ORANGE	1A920-6O
722500650	EASY RUN CAT. 6 - USA	1A920-6U
722500660	EASY RUN CAT. 6 - PURPLE	1A920-6V