

# HIGHLANDER®

Kategorivelger i forhold til brukervekt og aktivitetsnivå										
Aktivitet / Vekt	44 - 52 kg	53 - 59 kg	60 - 68 kg	69 - 77 kg	78 - 88 kg	89 - 100 kg	101 - 116 kg	117 - 130 kg	131 - 147 kg	148 - 166 kg
Lavt	1	1	2	3	4	5	6	7	8	9
Normalt	1	2	3	4	5	6	7	8	9	-
Høyt	2	3	4	5	6	7	8	9	-	-

Art.nr.	Modell	Lev.art.nr.
721101022	Highlander Foot Kit Cat 01 A Size 22 RU	FS3-00-01A22-RU
721101122	Highlander Foot Kit Cat 01 A Size 22 SL	FS3-00-01A22-SL
721101222	Highlander Foot Kit Cat 01 A Size 22 SR	FS3-00-01A22-SR
721101023	Highlander Foot Kit Cat 01 A Size 23 RU	FS3-00-01A23-RU
721101123	Highlander Foot Kit Cat 01 A Size 23 SL	FS3-00-01A23-SL
721101223	Highlander Foot Kit Cat 01 A Size 23 SR	FS3-00-01A23-SR
721101024	Highlander Foot Kit Cat 01 A Size 24 RU	FS3-00-01A24-RU
721101124	Highlander Foot Kit Cat 01 A Size 24 SL	FS3-00-01A24-SL
721101224	Highlander Foot Kit Cat 01 A Size 24 SR	FS3-00-01A24-SR
721101025	Highlander Foot Kit Cat 01 A Size 25 RU	FS3-00-01A25-RU
721101125	Highlander Foot Kit Cat 01 A Size 25 SL	FS3-00-01A25-SL
721101225	Highlander Foot Kit Cat 01 A Size 25 SR	FS3-00-01A25-SR
721101026	Highlander Foot Kit Cat 01 A Size 26 RU	FS3-00-01A26-RU
721101126	Highlander Foot Kit Cat 01 A Size 26 SL	FS3-00-01A26-SL
721101226	Highlander Foot Kit Cat 01 A Size 26 SR	FS3-00-01A26-SR
721101027	Highlander Foot Kit Cat 01 A Size 27 RU	FS3-00-01A27-RU
721101127	Highlander Foot Kit Cat 01 A Size 27 SL	FS3-00-01A27-SL
721101227	Highlander Foot Kit Cat 01 A Size 27 SR	FS3-00-01A27-SR
721101028	Highlander Foot Kit Cat 01 A Size 28 RU	FS3-00-01A28-RU
721101128	Highlander Foot Kit Cat 01 A Size 28 SL	FS3-00-01A28-SL
721101228	Highlander Foot Kit Cat 01 A Size 28 SR	FS3-00-01A28-SR
721101029	Highlander Foot Kit Cat 01 A Size 29 RU	FS3-00-01A29-RU
721101030	Highlander Foot Kit Cat 01 A Size 30 RU	FS3-00-01A30-RU
721101031	Highlander Foot Kit Cat 01 A Size 31 RU	FS3-00-01A31-RU
721102022	Highlander Foot Kit Cat 02 A Size 22 RU	FS3-00-02A22-RU
721102122	Highlander Foot Kit Cat 02 A Size 22 SL	FS3-00-02A22-SL
721102222	Highlander Foot Kit Cat 02 A Size 22 SR	FS3-00-02A22-SR
721102023	Highlander Foot Kit Cat 02 A Size 23 RU	FS3-00-02A23-RU
721102123	Highlander Foot Kit Cat 02 A Size 23 SL	FS3-00-02A23-SL
721102223	Highlander Foot Kit Cat 02 A Size 23 SR	FS3-00-02A23-SR
721102024	Highlander Foot Kit Cat 02 A Size 24 RU	FS3-00-02A24-RU
721102124	Highlander Foot Kit Cat 02 A Size 24 SL	FS3-00-02A24-SL
721102224	Highlander Foot Kit Cat 02 A Size 24 SR	FS3-00-02A24-SR
721102025	Highlander Foot Kit Cat 02 A Size 25 RU	FS3-00-02A25-RU
721102125	Highlander Foot Kit Cat 02 A Size 25 SL	FS3-00-02A25-SL
721102225	Highlander Foot Kit Cat 02 A Size 25 SR	FS3-00-02A25-SR

721102026	Highlander Foot Kit Cat 02 A Size 26 RU	FS3-00-02A26-RU
721102126	Highlander Foot Kit Cat 02 A Size 26 SL	FS3-00-02A26-SL
721102226	Highlander Foot Kit Cat 02 A Size 26 SR	FS3-00-02A26-SR
721102027	Highlander Foot Kit Cat 02 A Size 27 RU	FS3-00-02A27-RU
721102127	Highlander Foot Kit Cat 02 A Size 27 SL	FS3-00-02A27-SL
721102227	Highlander Foot Kit Cat 02 A Size 27 SR	FS3-00-02A27-SR
721102028	Highlander Foot Kit Cat 02 A Size 28 RU	FS3-00-02A28-RU
721102128	Highlander Foot Kit Cat 02 A Size 28 SL	FS3-00-02A28-SL
721102228	Highlander Foot Kit Cat 02 A Size 28 SR	FS3-00-02A28-SR
721102029	Highlander Foot Kit Cat 02 A Size 29 RU	FS3-00-02A29-RU
721102030	Highlander Foot Kit Cat 02 A Size 30 RU	FS3-00-02A30-RU
721102031	Highlander Foot Kit Cat 02 A Size 31 RU	FS3-00-02A31-RU
721103022	Highlander Foot Kit Cat 03 A Size 22 RU	FS3-00-03A22-RU
721103122	Highlander Foot Kit Cat 03 A Size 22 SL	FS3-00-03A22-SL
721103222	Highlander Foot Kit Cat 03 A Size 22 SR	FS3-00-03A22-SR
721103023	Highlander Foot Kit Cat 03 A Size 23 RU	FS3-00-03A23-RU
721103123	Highlander Foot Kit Cat 03 A Size 23 SL	FS3-00-03A23-SL
721103223	Highlander Foot Kit Cat 03 A Size 23 SR	FS3-00-03A23-SR
721103024	Highlander Foot Kit Cat 03 A Size 24 RU	FS3-00-03A24-RU
721103124	Highlander Foot Kit Cat 03 A Size 24 SL	FS3-00-03A24-SL
721103224	Highlander Foot Kit Cat 03 A Size 24 SR	FS3-00-03A24-SR
721103025	Highlander Foot Kit Cat 03 A Size 25 RU	FS3-00-03A25-RU
721103125	Highlander Foot Kit Cat 03 A Size 25 SL	FS3-00-03A25-SL
721103225	Highlander Foot Kit Cat 03 A Size 25 SR	FS3-00-03A25-SR
721103026	Highlander Foot Kit Cat 03 A Size 26 RU	FS3-00-03A26-RU
721103126	Highlander Foot Kit Cat 03 A Size 26 SL	FS3-00-03A26-SL
721103226	Highlander Foot Kit Cat 03 A Size 26 SR	FS3-00-03A26-SR
721103027	Highlander Foot Kit Cat 03 A Size 27 RU	FS3-00-03A27-RU
721103127	Highlander Foot Kit Cat 03 A Size 27 SL	FS3-00-03A27-SL
721103227	Highlander Foot Kit Cat 03 A Size 27 SR	FS3-00-03A27-SR
721103028	Highlander Foot Kit Cat 03 A Size 28 RU	FS3-00-03A28-RU
721103128	Highlander Foot Kit Cat 03 A Size 28 SL	FS3-00-03A28-SL
721103228	Highlander Foot Kit Cat 03 A Size 28 SR	FS3-00-03A28-SR
721103029	Highlander Foot Kit Cat 03 A Size 29 RU	FS3-00-03A29-RU
721103030	Highlander Foot Kit Cat 03 A Size 30 RU	FS3-00-03A30-RU
721103031	Highlander Foot Kit Cat 03 A Size 31 RU	FS3-00-03A31-RU
721104022	Highlander Foot Kit Cat 04 A Size 22 RU	FS3-00-04A22-RU
721104122	Highlander Foot Kit Cat 04 A Size 22 SL	FS3-00-04A22-SL
721104222	Highlander Foot Kit Cat 04 A Size 22 SR	FS3-00-04A22-SR
721104023	Highlander Foot Kit Cat 04 A Size 23 RU	FS3-00-04A23-RU
721104123	Highlander Foot Kit Cat 04 A Size 23 SL	FS3-00-04A23-SL
721104223	Highlander Foot Kit Cat 04 A Size 23 SR	FS3-00-04A23-SR
721104024	Highlander Foot Kit Cat 04 A Size 24 RU	FS3-00-04A24-RU
721104124	Highlander Foot Kit Cat 04 A Size 24 SL	FS3-00-04A24-SL
721104224	Highlander Foot Kit Cat 04 A Size 24 SR	FS3-00-04A24-SR

721104025	Highlander Foot Kit Cat 04 A Size 25 RU	FS3-00-04A25-RU
721104125	Highlander Foot Kit Cat 04 A Size 25 SL	FS3-00-04A25-SL
721104225	Highlander Foot Kit Cat 04 A Size 25 SR	FS3-00-04A25-SR
721104026	Highlander Foot Kit Cat 04 A Size 26 RU	FS3-00-04A26-RU
721104126	Highlander Foot Kit Cat 04 A Size 26 SL	FS3-00-04A26-SL
721104226	Highlander Foot Kit Cat 04 A Size 26 SR	FS3-00-04A26-SR
721104027	Highlander Foot Kit Cat 04 A Size 27 RU	FS3-00-04A27-RU
721104127	Highlander Foot Kit Cat 04 A Size 27 SL	FS3-00-04A27-SL
721104227	Highlander Foot Kit Cat 04 A Size 27 SR	FS3-00-04A27-SR
721104028	Highlander Foot Kit Cat 04 A Size 28 RU	FS3-00-04A28-RU
721104128	Highlander Foot Kit Cat 04 A Size 28 SL	FS3-00-04A28-SL
721104228	Highlander Foot Kit Cat 04 A Size 28 SR	FS3-00-04A28-SR
721104029	Highlander Foot Kit Cat 04 A Size 29 RU	FS3-00-04A29-RU
721104030	Highlander Foot Kit Cat 04 A Size 30 RU	FS3-00-04A30-RU
721104031	Highlander Foot Kit Cat 04 A Size 31 RU	FS3-00-04A31-RU
721105022	Highlander Foot Kit Cat 05 A Size 22 RU	FS3-00-05A22-RU
721105122	Highlander Foot Kit Cat 05 A Size 22 SL	FS3-00-05A22-SL
721105222	Highlander Foot Kit Cat 05 A Size 22 SR	FS3-00-05A22-SR
721105023	Highlander Foot Kit Cat 05 A Size 23 RU	FS3-00-05A23-RU
721105123	Highlander Foot Kit Cat 05 A Size 23 SL	FS3-00-05A23-SL
721105223	Highlander Foot Kit Cat 05 A Size 23 SR	FS3-00-05A23-SR
721105024	Highlander Foot Kit Cat 05 A Size 24 RU	FS3-00-05A24-RU
721105124	Highlander Foot Kit Cat 05 A Size 24 SL	FS3-00-05A24-SL
721105224	Highlander Foot Kit Cat 05 A Size 24 SR	FS3-00-05A24-SR
721105025	Highlander Foot Kit Cat 05 A Size 25 RU	FS3-00-05A25-RU
721105125	Highlander Foot Kit Cat 05 A Size 25 SL	FS3-00-05A25-SL
721105225	Highlander Foot Kit Cat 05 A Size 25 SR	FS3-00-05A25-SR
721105026	Highlander Foot Kit Cat 05 A Size 26 RU	FS3-00-05A26-RU
721105126	Highlander Foot Kit Cat 05 A Size 26 SL	FS3-00-05A26-SL
721105226	Highlander Foot Kit Cat 05 A Size 26 SR	FS3-00-05A26-SR
721105027	Highlander Foot Kit Cat 05 A Size 27 RU	FS3-00-05A27-RU
721105127	Highlander Foot Kit Cat 05 A Size 27 SL	FS3-00-05A27-SL
721105227	Highlander Foot Kit Cat 05 A Size 27 SR	FS3-00-05A27-SR
721105028	Highlander Foot Kit Cat 05 A Size 28 RU	FS3-00-05A28-RU
721105128	Highlander Foot Kit Cat 05 A Size 28 SL	FS3-00-05A28-SL
721105228	Highlander Foot Kit Cat 05 A Size 28 SR	FS3-00-05A28-SR
721105029	Highlander Foot Kit Cat 05 A Size 29 RU	FS3-00-05A29-RU
721105030	Highlander Foot Kit Cat 05 A Size 30 RU	FS3-00-05A30-RU
721105031	Highlander Foot Kit Cat 05 A Size 31 RU	FS3-00-05A31-RU
721106022	Highlander Foot Kit Cat 06 A Size 22 RU	FS3-00-06A22-RU
721106122	Highlander Foot Kit Cat 06 A Size 22 SL	FS3-00-06A22-SL
721106222	Highlander Foot Kit Cat 06 A Size 22 SR	FS3-00-06A22-SR
721106023	Highlander Foot Kit Cat 06 A Size 23 RU	FS3-00-06A23-RU
721106123	Highlander Foot Kit Cat 06 A Size 23 SL	FS3-00-06A23-SL
721106223	Highlander Foot Kit Cat 06 A Size 23 SR	FS3-00-06A23-SR

721106024	Highlander Foot Kit Cat 06 A Size 24 RU	FS3-00-06A24-RU
721106124	Highlander Foot Kit Cat 06 A Size 24 SL	FS3-00-06A24-SL
721106224	Highlander Foot Kit Cat 06 A Size 24 SR	FS3-00-06A24-SR
721106025	Highlander Foot Kit Cat 06 A Size 25 RU	FS3-00-06A25-RU
721106125	Highlander Foot Kit Cat 06 A Size 25 SL	FS3-00-06A25-SL
721106225	Highlander Foot Kit Cat 06 A Size 25 SR	FS3-00-06A25-SR
721106026	Highlander Foot Kit Cat 06 A Size 26 RU	FS3-00-06A26-RU
721106126	Highlander Foot Kit Cat 06 A Size 26 SL	FS3-00-06A26-SL
721106226	Highlander Foot Kit Cat 06 A Size 26 SR	FS3-00-06A26-SR
721106027	Highlander Foot Kit Cat 06 A Size 27 RU	FS3-00-06A27-RU
721106127	Highlander Foot Kit Cat 06 A Size 27 SL	FS3-00-06A27-SL
721106227	Highlander Foot Kit Cat 06 A Size 27 SR	FS3-00-06A27-SR
721106028	Highlander Foot Kit Cat 06 A Size 28 RU	FS3-00-06A28-RU
721106128	Highlander Foot Kit Cat 06 A Size 28 SL	FS3-00-06A28-SL
721106228	Highlander Foot Kit Cat 06 A Size 28 SR	FS3-00-06A28-SR
721106029	Highlander Foot Kit Cat 06 A Size 29 RU	FS3-00-06A29-RU
721106030	Highlander Foot Kit Cat 06 A Size 30 RU	FS3-00-06A30-RU
721106031	Highlander Foot Kit Cat 06 A Size 31 RU	FS3-00-06A31-RU
721107022	Highlander Foot Kit Cat 07 A Size 22 RU	FS3-00-07A22-RU
721107122	Highlander Foot Kit Cat 07 A Size 22 SL	FS3-00-07A22-SL
721107222	Highlander Foot Kit Cat 07 A Size 22 SR	FS3-00-07A22-SR
721107023	Highlander Foot Kit Cat 07 A Size 23 RU	FS3-00-07A23-RU
721107123	Highlander Foot Kit Cat 07 A Size 23 SL	FS3-00-07A23-SL
721107223	Highlander Foot Kit Cat 07 A Size 23 SR	FS3-00-07A23-SR
721107024	Highlander Foot Kit Cat 07 A Size 24 RU	FS3-00-07A24-RU
721107124	Highlander Foot Kit Cat 07 A Size 24 SL	FS3-00-07A24-SL
721107224	Highlander Foot Kit Cat 07 A Size 24 SR	FS3-00-07A24-SR
721107025	Highlander Foot Kit Cat 07 A Size 25 RU	FS3-00-07A25-RU
721107125	Highlander Foot Kit Cat 07 A Size 25 SL	FS3-00-07A25-SL
721107225	Highlander Foot Kit Cat 07 A Size 25 SR	FS3-00-07A25-SR
721107026	Highlander Foot Kit Cat 07 A Size 26 RU	FS3-00-07A26-RU
721107126	Highlander Foot Kit Cat 07 A Size 26 SL	FS3-00-07A26-SL
721107226	Highlander Foot Kit Cat 07 A Size 26 SR	FS3-00-07A26-SR
721107027	Highlander Foot Kit Cat 07 A Size 27 RU	FS3-00-07A27-RU
721107127	Highlander Foot Kit Cat 07 A Size 27 SL	FS3-00-07A27-SL
721107227	Highlander Foot Kit Cat 07 A Size 27 SR	FS3-00-07A27-SR
721107028	Highlander Foot Kit Cat 07 A Size 28 RU	FS3-00-07A28-RU
721107128	Highlander Foot Kit Cat 07 A Size 28 SL	FS3-00-07A28-SL
721107228	Highlander Foot Kit Cat 07 A Size 28 SR	FS3-00-07A28-SR
721107029	Highlander Foot Kit Cat 07 A Size 29 RU	FS3-00-07A29-RU
721107030	Highlander Foot Kit Cat 07 A Size 30 RU	FS3-00-07A30-RU
721107031	Highlander Foot Kit Cat 07 A Size 31 RU	FS3-00-07A31-RU
721108022	Highlander Foot Kit Cat 08 A Size 22 RU	FS3-00-08A22-RU
721108122	Highlander Foot Kit Cat 08 A Size 22 SL	FS3-00-08A22-SL
721108222	Highlander Foot Kit Cat 08 A Size 22 SR	FS3-00-08A22-SR

721108023	Highlander Foot Kit Cat 08 A Size 23	RU	FS3-00-08A23-RU
721108123	Highlander Foot Kit Cat 08 A Size 23	SL	FS3-00-08A23-SL
721108223	Highlander Foot Kit Cat 08 A Size 23	SR	FS3-00-08A23-SR
721108024	Highlander Foot Kit Cat 08 A Size 24	RU	FS3-00-08A24-RU
721108124	Highlander Foot Kit Cat 08 A Size 24	SL	FS3-00-08A24-SL
721108224	Highlander Foot Kit Cat 08 A Size 24	SR	FS3-00-08A24-SR
721108025	Highlander Foot Kit Cat 08 A Size 25	RU	FS3-00-08A25-RU
721108125	Highlander Foot Kit Cat 08 A Size 25	SL	FS3-00-08A25-SL
721108225	Highlander Foot Kit Cat 08 A Size 25	SR	FS3-00-08A25-SR
721108026	Highlander Foot Kit Cat 08 A Size 26	RU	FS3-00-08A26-RU
721108126	Highlander Foot Kit Cat 08 A Size 26	SL	FS3-00-08A26-SL
721108226	Highlander Foot Kit Cat 08 A Size 26	SR	FS3-00-08A26-SR
721108027	Highlander Foot Kit Cat 08 A Size 27	RU	FS3-00-08A27-RU
721108127	Highlander Foot Kit Cat 08 A Size 27	SL	FS3-00-08A27-SL
721108227	Highlander Foot Kit Cat 08 A Size 27	SR	FS3-00-08A27-SR
721108028	Highlander Foot Kit Cat 08 A Size 28	RU	FS3-00-08A28-RU
721108128	Highlander Foot Kit Cat 08 A Size 28	SL	FS3-00-08A28-SL
721108228	Highlander Foot Kit Cat 08 A Size 28	SR	FS3-00-08A28-SR
721108029	Highlander Foot Kit Cat 08 A Size 29	RU	FS3-00-08A29-RU
721108030	Highlander Foot Kit Cat 08 A Size 30	RU	FS3-00-08A30-RU
721108031	Highlander Foot Kit Cat 08 A Size 31	RU	FS3-00-08A31-RU
721109022	Highlander Foot Kit Cat 09 A Size 22	RU	FS3-00-09A22-RU
721109122	Highlander Foot Kit Cat 09 A Size 22	SL	FS3-00-09A22-SL
721109222	Highlander Foot Kit Cat 09 A Size 22	SR	FS3-00-09A22-SR
721109023	Highlander Foot Kit Cat 09 A Size 23	RU	FS3-00-09A23-RU
721109123	Highlander Foot Kit Cat 09 A Size 23	SL	FS3-00-09A23-SL
721109223	Highlander Foot Kit Cat 09 A Size 23	SR	FS3-00-09A23-SR
721109024	Highlander Foot Kit Cat 09 A Size 24	RU	FS3-00-09A24-RU
721109124	Highlander Foot Kit Cat 09 A Size 24	SL	FS3-00-09A24-SL
721109224	Highlander Foot Kit Cat 09 A Size 24	SR	FS3-00-09A24-SR
721109025	Highlander Foot Kit Cat 09 A Size 25	RU	FS3-00-09A25-RU
721109125	Highlander Foot Kit Cat 09 A Size 25	SL	FS3-00-09A25-SL
721109225	Highlander Foot Kit Cat 09 A Size 25	SR	FS3-00-09A25-SR
721109026	Highlander Foot Kit Cat 09 A Size 26	RU	FS3-00-09A26-RU
721109126	Highlander Foot Kit Cat 09 A Size 26	SL	FS3-00-09A26-SL
721109226	Highlander Foot Kit Cat 09 A Size 26	SR	FS3-00-09A26-SR
721109027	Highlander Foot Kit Cat 09 A Size 27	RU	FS3-00-09A27-RU
721109127	Highlander Foot Kit Cat 09 A Size 27	SL	FS3-00-09A27-SL
721109227	Highlander Foot Kit Cat 09 A Size 27	SR	FS3-00-09A27-SR
721109028	Highlander Foot Kit Cat 09 A Size 28	RU	FS3-00-09A28-RU
721109128	Highlander Foot Kit Cat 09 A Size 28	SL	FS3-00-09A28-SL
721109228	Highlander Foot Kit Cat 09 A Size 28	SR	FS3-00-09A28-SR
721109029	Highlander Foot Kit Cat 09 A Size 29	RU	FS3-00-09A29-RU
721109030	Highlander Foot Kit Cat 09 A Size 30	RU	FS3-00-09A30-RU
721109031	Highlander Foot Kit Cat 09 A Size 31	RU	FS3-00-09A31-RU

---