

KINNEX 2.0

Kategorivelger i forhold til aktivitetsnivå, brukervekt og størrelse								
Vekt / Aktivitetsnivå	44-52 kg	53-59 kg	60-68 kg	69-77 kg	78-88 kg	89-100 kg	101-116 kg	117-132 kg
Lavt	1	1	2	3	4	5	6	7
Moderat	1	2	3	4	5	6	7	8

Art.nr.	Modell	Lev.art.nr.
722270124	Kinnex 2.0 Foot/Ankle, Cat.1 Size 24cm	F14-N2-01A24-RU
722271124	Kinnex 2.0 Foot/Ankle, Cat. 1 Size 24cm Sandal Left	F14-N2-01A24-SL
722272124	Kinnex 2.0 Foot/Ankle, Cat. 1 Size 24cm Sandal Right	F14-N2-01A24-SR
722270125	Kinnex 2.0 Foot/Ankle, Cat. 1 Size 25cm	F14-N2-01A25-RU
722271125	Kinnex 2.0 Foot/Ankle, Cat. 1 Size 25cm Sandal Left	F14-N2-01A25-SL
722272125	Kinnex 2.0 Foot/Ankle, Cat. 1 Size 25cm Sandal Right	F14-N2-01A25-SR
722270126	Kinnex 2.0 Foot/Ankle, Cat. 1 Size 26cm	F14-N2-01A26-RU
722271126	Kinnex 2.0 Foot/Ankle, Cat. 1 Size 26cm Sandal Left	F14-N2-01A26-SL
722272126	Kinnex 2.0 Foot/Ankle, Cat. 1 Size 26cm Sandal Right	F14-N2-01A26-SR
722270127	Kinnex 2.0 Foot/Ankle, Cat. 1 Size 27cm	F14-N2-01A27-RU
722271127	Kinnex 2.0 Foot/Ankle, Cat. 1 Size 27cm Sandal Left	F14-N2-01A27-SL
722272127	Kinnex 2.0 Foot/Ankle, Cat. 1 Size 27cm Sandal Right	F14-N2-01A27-SR
722270128	Kinnex 2.0 Foot/Ankle, Cat. 1 Size 28cm	F14-N2-01A28-RU
722271128	Kinnex 2.0 Foot/Ankle, Cat. 1 Size 28cm Sandal Left	F14-N2-01A28-SL
722272128	Kinnex 2.0 Foot/Ankle, Cat. 1 Size 28cm Sandal Right	F14-N2-01A28-SR
722270129	Kinnex 2.0 Foot/Ankle, Cat. 1 Size 29cm	F14-N2-01A29-RU
722270130	Kinnex 2.0 Foot/Ankle, Cat. 1 Size 30cm	F14-N2-01A30-RU
722270224	Kinnex 2.0 Foot/Ankle, Cat. 2 Size 24cm	F14-N2-02A24-RU
722271224	Kinnex 2.0 Foot/Ankle, Cat. 2 Size 24cm Sandal Left	F14-N2-02A24-SL
722272224	Kinnex 2.0 Foot/Ankle, Cat. 2 Size 24cm Sandal Right	F14-N2-02A24-SR
722270225	Kinnex 2.0 Foot/Ankle, Cat. 2 Size 25cm	F14-N2-02A25-RU
722271225	Kinnex 2.0 Foot/Ankle, Cat. 2 Size 25cm Sandal Left	F14-N2-02A25-SL
722272225	Kinnex 2.0 Foot/Ankle, Cat. 2 Size 25cm Sandal Right	F14-N2-02A25-SR
722270226	Kinnex 2.0 Foot/Ankle, Cat. 2 Size 26cm	F14-N2-02A26-RU
722271226	Kinnex 2.0 Foot/Ankle, Cat. 2 Size 26cm Sandal Left	F14-N2-02A26-SL
722272226	Kinnex 2.0 Foot/Ankle, Cat. 2 Size 26cm Sandal Right	F14-N2-02A26-SR
722270227	Kinnex 2.0 Foot/Ankle, Cat. 2 Size 27cm	F14-N2-02A27-RU
722271227	Kinnex 2.0 Foot/Ankle, Cat. 2 Size 27cm Sandal Left	F14-N2-02A27-SL
722272227	Kinnex 2.0 Foot/Ankle, Cat. 2 Size 27cm Sandal Right	F14-N2-02A27-SR
722270228	Kinnex 2.0 Foot/Ankle, Cat. 2 Size 28cm	F14-N2-02A28-RU
722271228	Kinnex 2.0 Foot/Ankle, Cat. 2 Size 28cm Sandal Left	F14-N2-02A28-SL
722272228	Kinnex 2.0 Foot/Ankle, Cat. 2 Size 28cm Sandal Right	F14-N2-02A28-SR
722270229	Kinnex 2.0 Foot/Ankle, Cat. 2 Size 29cm	F14-N2-02A29-RU
722270230	Kinnex 2.0 Foot/Ankle, Cat. 2 Size 30cm	F14-N2-02A30-RU
722270324	Kinnex 2.0 Foot/Ankle, Cat. 3 Size 24cm	F14-N2-03A24-RU
722271324	Kinnex 2.0 Foot/Ankle, Cat. 3 Size 24cm Sandal Left	F14-N2-03A24-SL
722272324	Kinnex 2.0 Foot/Ankle, Cat. 3 Size 24cm Sandal Right	F14-N2-03A24-SR
722270325	Kinnex 2.0 Foot/Ankle, Cat. 3 Size 25cm	F14-N2-03A25-RU
722271325	Kinnex 2.0 Foot/Ankle, Cat. 3 Size 25cm Sandal Left	F14-N2-03A25-SL
722272325	Kinnex 2.0 Foot/Ankle, Cat. 3 Size 25cm Sandal Right	F14-N2-03A25-SR
722270326	Kinnex 2.0 Foot/Ankle, Cat. 3 Size 26cm	F14-N2-03A26-RU
722271326	Kinnex 2.0 Foot/Ankle, Cat. 3 Size 26cm Sandal Left	F14-N2-03A26-SL
722272326	Kinnex 2.0 Foot/Ankle, Cat. 3 Size 26cm Sandal Right	F14-N2-03A26-SR

722270327	Kinnex 2.0 Foot/Ankle, Cat. 3 Size 27cm	F14-N2-03A27-RU
722271327	Kinnex 2.0 Foot/Ankle, Cat. 3 Size 27cm Sandal Left	F14-N2-03A27-SL
722272327	Kinnex 2.0 Foot/Ankle, Cat. 3 Size 27cm Sandal Right	F14-N2-03A27-SR
722270328	Kinnex 2.0 Foot/Ankle, Cat. 3 Size 28cm	F14-N2-03A28-RU
722271328	Kinnex 2.0 Foot/Ankle, Cat. 3 Size 28cm Sandal Left	F14-N2-03A28-SL
722272328	Kinnex 2.0 Foot/Ankle, Cat. 3 Size 28cm Sandal Right	F14-N2-03A28-SR
722270329	Kinnex 2.0 Foot/Ankle, Cat. 3 Size 29cm	F14-N2-03A29-RU
722270330	Kinnex 2.0 Foot/Ankle, Cat. 3 Size 30cm	F14-N2-03A30-RU
722270424	Kinnex 2.0 Foot/Ankle, Cat. 4 Size 24cm	F14-N2-04A24-RU
722271424	Kinnex 2.0 Foot/Ankle, Cat. 4 Size 24cm Sandal Left	F14-N2-04A24-SL
722272424	Kinnex 2.0 Foot/Ankle, Cat. 4 Size 24cm Sandal Right	F14-N2-04A24-SR
722270425	Kinnex 2.0 Foot/Ankle, Cat. 4 Size 25cm	F14-N2-04A25-RU
722271425	Kinnex 2.0 Foot/Ankle, Cat. 4 Size 25cm Sandal Left	F14-N2-04A25-SL
722272425	Kinnex 2.0 Foot/Ankle, Cat. 4 Size 25cm Sandal Right	F14-N2-04A25-SR
722270426	Kinnex 2.0 Foot/Ankle, Cat. 4 Size 26cm	F14-N2-04A26-RU
722271426	Kinnex 2.0 Foot/Ankle, Cat. 4 Size 26cm Sandal Left	F14-N2-04A26-SL
722272426	Kinnex 2.0 Foot/Ankle, Cat. 4 Size 26cm Sandal Right	F14-N2-04A26-SR
722270427	Kinnex 2.0 Foot/Ankle, Cat. 4 Size 27cm	F14-N2-04A27-RU
722271427	Kinnex 2.0 Foot/Ankle, Cat. 4 Size 27cm Sandal Left	F14-N2-04A27-SL
722272427	Kinnex 2.0 Foot/Ankle, Cat. 4 Size 27cm Sandal Right	F14-N2-04A27-SR
722270428	Kinnex 2.0 Foot/Ankle, Cat. 4 Size 28cm	F14-N2-04A28-RU
722271428	Kinnex 2.0 Foot/Ankle, Cat. 4 Size 28cm Sandal Left	F14-N2-04A28-SL
722272428	Kinnex 2.0 Foot/Ankle, Cat. 4 Size 28cm Sandal Right	F14-N2-04A28-SR
722270429	Kinnex 2.0 Foot/Ankle, Cat. 4 Size 29cm	F14-N2-04A29-RU
722270430	Kinnex 2.0 Foot/Ankle, Cat. 4 Size 30cm	F14-N2-04A30-RU
722270524	Kinnex 2.0 Foot/Ankle, Cat. 5 Size 24cm	F14-N2-05A24-RU
722271524	Kinnex 2.0 Foot/Ankle, Cat. 5 Size 24cm Sandal Left	F14-N2-05A24-SL
722272524	Kinnex 2.0 Foot/Ankle, Cat. 5 Size 24cm Sandal Right	F14-N2-05A24-SR
722270525	Kinnex 2.0 Foot/Ankle, Cat. 5 Size 25cm	F14-N2-05A25-RU
722271525	Kinnex 2.0 Foot/Ankle, Cat. 5 Size 25cm Sandal Left	F14-N2-05A25-SL
722272525	Kinnex 2.0 Foot/Ankle, Cat. 5 Size 25cm Sandal Right	F14-N2-05A25-SR
722270526	Kinnex 2.0 Foot/Ankle, Cat. 5 Size 26cm	F14-N2-05A26-RU
722271526	Kinnex 2.0 Foot/Ankle, Cat. 5 Size 26cm Sandal Left	F14-N2-05A26-SL
722272526	Kinnex 2.0 Foot/Ankle, Cat. 5 Size 26cm Sandal Right	F14-N2-05A26-SR
722270527	Kinnex 2.0 Foot/Ankle, Cat. 5 Size 27cm	F14-N2-05A27-RU
722271527	Kinnex 2.0 Foot/Ankle, Cat. 5 Size 27cm Sandal Left	F14-N2-05A27-SL
722272527	Kinnex 2.0 Foot/Ankle, Cat. 5 Size 27cm Sandal Right	F14-N2-05A27-SR
722270528	Kinnex 2.0 Foot/Ankle, Cat. 5 Size 28cm	F14-N2-05A28-RU
722271528	Kinnex 2.0 Foot/Ankle, Cat. 5 Size 28cm Sandal Left	F14-N2-05A28-SL
722272528	Kinnex 2.0 Foot/Ankle, Cat. 5 Size 28cm Sandal Right	F14-N2-05A28-SR
722270529	Kinnex 2.0 Foot/Ankle, Cat. 5 Size 29cm	F14-N2-05A29-RU
722270530	Kinnex 2.0 Foot/Ankle, Cat. 5 Size 30cm	F14-N2-05A30-RU
722270624	Kinnex 2.0 Foot/Ankle, Cat. 6 Size 24cm	F14-N2-06A24-RU
722271624	Kinnex 2.0 Foot/Ankle, Cat. 6 Size 24cm Sandal Left	F14-N2-06A24-SL
722272624	Kinnex 2.0 Foot/Ankle, Cat. 6 Size 24cm Sandal Right	F14-N2-06A24-SR
722270625	Kinnex 2.0 Foot/Ankle, Cat. 6 Size 25cm	F14-N2-06A25-RU
722271625	Kinnex 2.0 Foot/Ankle, Cat. 6 Size 25cm Sandal Left	F14-N2-06A25-SL
722272625	Kinnex 2.0 Foot/Ankle, Cat. 6 Size 25cm Sandal Right	F14-N2-06A25-SR
722270626	Kinnex 2.0 Foot/Ankle, Cat. 6 Size 26cm	F14-N2-06A26-RU
722271626	Kinnex 2.0 Foot/Ankle, Cat. 6 Size 26cm Sandal Left	F14-N2-06A26-SL
722272626	Kinnex 2.0 Foot/Ankle, Cat. 6 Size 26cm Sandal Right	F14-N2-06A26-SR

722270627	Kinnex 2.0 Foot/Ankle, Cat. 6 Size 27cm	F14-N2-06A27-RU
722271627	Kinnex 2.0 Foot/Ankle, Cat. 6 Size 27cm Sandal Left	F14-N2-06A27-SL
722272627	Kinnex 2.0 Foot/Ankle, Cat. 6 Size 27cm Sandal Right	F14-N2-06A27-SR
722270628	Kinnex 2.0 Foot/Ankle, Cat. 6 Size 28cm	F14-N2-06A28-RU
722271628	Kinnex 2.0 Foot/Ankle, Cat. 6 Size 28cm Sandal Left	F14-N2-06A28-SL
722272628	Kinnex 2.0 Foot/Ankle, Cat. 6 Size 28cm Sandal Right	F14-N2-06A28-SR
722270629	Kinnex 2.0 Foot/Ankle, Cat. 6 Size 29cm	F14-N2-06A29-RU
722270630	Kinnex 2.0 Foot/Ankle, Cat. 6 Size 30cm	F14-N2-06A30-RU
722270724	Kinnex 2.0 Foot/Ankle, Cat. 7 Size 24cm	F14-N2-07A24-RU
722271724	Kinnex 2.0 Foot/Ankle, Cat. 7 Size 24cm Sandal Left	F14-N2-07A24-SL
722272724	Kinnex 2.0 Foot/Ankle, Cat. 7 Size 24cm Sandal Right	F14-N2-07A24-SR
722270725	Kinnex 2.0 Foot/Ankle, Cat. 7 Size 25cm	F14-N2-07A25-RU
722271725	Kinnex 2.0 Foot/Ankle, Cat. 7 Size 25cm Sandal Left	F14-N2-07A25-SL
722272725	Kinnex 2.0 Foot/Ankle, Cat. 7 Size 25cm Sandal Right	F14-N2-07A25-SR
722270726	Kinnex 2.0 Foot/Ankle, Cat. 7 Size 26cm	F14-N2-07A26-RU
722271726	Kinnex 2.0 Foot/Ankle, Cat. 7 Size 26cm Sandal Left	F14-N2-07A26-SL
722272726	Kinnex 2.0 Foot/Ankle, Cat. 7 Size 26cm Sandal Right	F14-N2-07A26-SR
722270727	Kinnex 2.0 Foot/Ankle, Cat. 7 Size 27cm	F14-N2-07A27-RU
722271727	Kinnex 2.0 Foot/Ankle, Cat. 7 Size 27cm Sandal Left	F14-N2-07A27-SL
722272727	Kinnex 2.0 Foot/Ankle, Cat. 7 Size 27cm Sandal Right	F14-N2-07A27-SR
722270728	Kinnex 2.0 Foot/Ankle, Cat. 7 Size 28cm	F14-N2-07A28-RU
722271728	Kinnex 2.0 Foot/Ankle, Cat. 7 Size 28cm Sandal Left	F14-N2-07A28-SL
722272728	Kinnex 2.0 Foot/Ankle, Cat. 7 Size 28cm Sandal Right	F14-N2-07A28-SR
722270729	Kinnex 2.0 Foot/Ankle, Cat. 7 Size 29cm	F14-N2-07A29-RU
722270730	Kinnex 2.0 Foot/Ankle, Cat. 7 Size 30cm	F14-N2-07A30-RU