

SPECIALTY

Highlander[®] MAX

Proven technology for dependable performance

- Offers an ideal blend of stability and action, often preferred by transfemoral amputees
- Split toe provides excellent inversion/eversion
- Custom Highlander MAX available for users weighing up to 226 kg, sold with optional heavy duty 34 mm component kit
- 36-month warranty



Split Toe




Sandal Toe

Highlander MAX (FS3-H5)

Highlander® MAX



SPECIALTY



Note: Build height measurements represent a fully assembled unit.

Product ID: FS3-H5

Sizes: 26-31 cm

Sandal toe available for sizes 22-28 cm*

User Weight Rating: Rated to 226 kg with 34 mm pyramid**

Heel Height: 10 mm

Connector: Pyramid

Warranty: Foot module 36 months, Cosmetic foot shell 6 months

* When ordering the Sandal Toe option, please contact Customer Service for lead time.

**Custom orders: 7-10 day lead time, no trial period.

Accessories

- 34 mm Component Kit: (KIT-00-202), rated to 226 kg; kit includes pylon, two female tube clamps, and a 4-ear connector
 - Tube Clamp, 34 mm: (ACC-00-20201-00)
 - 4-Ear Connector: (ACC-00-20202-00)
 - Pylon, 34 mm: (ACC-00-20203-00)
- Foot Shell Removal Tool: (ACC-00-103)
- Spectra Sock: (ACC-00-101xx)
30 cm, 35 cm, 40 cm
- Stiffening Bumpers:
 - 22-25 cm (KIT-00-11303)
 - 26-28 cm (KIT-00-11304)
 - 29-31 cm (KIT-00-11305)

Note: Spectra socks and heel bumpers are included in the delivery.

Build Height & Weight

Size (cm)	Height (mm)	Weight (g)
26-28	121	590
29-31	130	750

Product weight includes foot shell. Measurements are approximate and will vary due to category selection and manufacturing variability.

Stiffness Category Selection

Impact Level			Weight (kg)
Low	Moderate	High	
		10	131-147
	10	11	148-166
10			167-190
11			191-218
12			219-226

Complete Prosthetic Foot (incl. foot shell) Part Number

FS3-H5-A-M-3

Category 10-12 Size 22-31 Side L-left R-right Skin Tone L-light M-medium D-dark

Foot Module Part Number

FS3-H5-A-

Category 10-12 Size 22-31 Toe Configuration RU-regular toe, universal SR-sandal toe, right SL-sandal toe, left

Foot Shell Part Number

FTC--4-

Generation 2-3 Profile M-mid H-high Cap 0=yes 1=no Size 22-31 Skin Tone L-light M-medium D-dark Toe R-regular S-sandal Side R-right L-left

See the Cosmetic Foot Shell page for more detailed ordering instructions.