


Realistic Prosthetics Ltd.

Dimension Chart for **Passive** Cosmetic Gloves.

GLOVE	Hand / Gender	A Across Knuckles	B Palm Circumference	C Wrist Circumference	D Middle Finger	Length to Thumb Tip	Remarks
RPL 501	Left Male	3 ½" 89 mm	220 mm	185 mm	78 mm	575 mm	
RPL 502	Right Male	3 ½" 89 mm	220 mm	185 mm	78 mm	575 mm	
RPL 503	Left Male	3 ¼" 82.5 mm	186 mm	146 mm	71 mm	520 mm	
RPL 504	Right Male	3 ¼" 82.5 mm	186 mm	146 mm	71 mm	520 mm	
RPL 521	Left Male	3 ¼" 82.5 mm	197 mm	168 mm	78 mm	550 mm	Long Fingers
RPL 522	Right Male	3 ¼" 82.5 mm	197 mm	168 mm	78 mm	550 mm	Long Fingers
RPL 505	Left Female	3 " 76 mm	182 mm	150 mm	70 mm	510 mm	
RPL 506	Right Female	3 " 76 mm	182 mm	150 mm	70 mm	510 mm	
RPL 505M	Left Male	3 " 76 mm	182 mm	150 mm	69 mm	490 mm	
RPL 506M	Right Male	3 " 76 mm	182 mm	150 mm	69 mm	490 mm	
RPL 507	Left Female	2 ¾" 70 mm	168 mm	135 mm	70 mm	500 mm	
RPL 508	Right Female	2 ¾" 70 mm	168 mm	135 mm	70 mm	500 mm	
RPL 511	Left Child	2 ½" 64 mm	152 mm	128 mm	56 mm	420 mm	
RPL 512	Right Child	2 ½" 64 mm	152 mm	128 mm	56 mm	420 mm	
RPL 513	Left Child	2 ¼" 57 mm	132 mm	126 mm	48 mm	340 mm	
RPL 514	Right Child	2 ¼" 57 mm	132 mm	126 mm	48 mm	340 mm	



Our cosmetic gloves are compatible with UK, European and other manufacturers hand systems
All our gloves are CE marked

Realistic Prosthetics Ltd.

Unit 19 Platts Eyot, Lower Sunbury Road, Hampton, Middlesex TW12 2HF U.K.

Phone: 020 8783 0099

Fax: 020 8941 4077

International +44 (0)

Web: www.realisticprosthetics.com

E-mail: info@realisticprosthetics.com