



REALISTIC PROSTHETICS

Myo/Mech. Realskin Gloves

Lifelike
appearance
RealSkin™



Choice of realistic
skin tones

Compatible with
UK and other hand
systems





Hand finished nails
give the gloves their truly
lifelike appearance



Cosmetic Glove Specialists

Extra length fitting
above elbow

Dimension Chart for Myo/Mech. Cosmetic Gloves



GLOVE	Otto Bock Size*	Hand	Across Knuckles	A Palm Circumference	B Wrist Circumference	C Middle Finger	D Length to Thumb tip
RPL 401	8¼	Left	89mm	219mm	176mm	76mm	580mm
RPL 402	8¼	Right	89mm	219mm	176mm	76mm	580mm
RPL 403	7¾	Left	82.5mm	205mm	168mm	71mm	570mm
RPL 404	7¾	Right	82.5mm	205mm	168mm	71mm	570mm
RPL 405	7¼	Left	76mm	194mm	168mm	71mm	540mm
RPL 406	7¼	Right	76mm	194mm	168mm	71mm	540mm
RPL 407	6¾	Left	70mm	174mm	150mm	67mm	440mm
RPL 408	6¾	Right	70mm	174mm	150mm	67mm	440mm

*Otto Bock is the Trade Mark of Otto Bock Healthcare Plc.

Realistic Prosthetics Ltd
 Unit 19 Platts Eyot, Lower Sunbury Road,
 Hampton, Middlesex TW12 2HF UK
 Phone +44 (0) 20 8783 0099
 Fax +44 (0) 20 8941 4077
 E-mail info@realisticprosthetics.com
realisticprosthetics.com

These cosmetic gloves are compatible with UK, European and other manufacturers hand systems

CE marked product

For sales, please contact us directly, or our distribution Agent in your Country.