



REALISTIC PROSTHETICS

RealSkin Extra Length Gloves

Lifelike real skin
appearance



Choice of realistic
skin tones

Flexible & rigid
options available




Hand finished nails
give the gloves their truly
lifelike appearance



Cosmetic Glove Specialists

Extra length fitting
above elbow

Dimension Chart for Passive Cosmetic Gloves



GLOVE	Hand / Gender	Across Knuckles		A	B	C	D
		inch	mm	Palm Circumference	Wrist Circumference	Middle Finger	Length to Thumb tip
RPL 501	Left Male	3½"	89mm	220mm	185mm	78mm	575mm
RPL 502	Right Male	3½"	89mm	220mm	185mm	78mm	575mm
RPL 503	Left Male	3¼"	82.5mm	186mm	146mm	71mm	520mm
RPL 504	Right Male	3¼"	82.5mm	186mm	146mm	71mm	520mm
RPL 521	Left Male LF	3¼"	82.5mm	197mm	168mm	78mm	550mm
RPL 522	Right Male LF	3¼"	82.5mm	197mm	168mm	78mm	550mm
RPL 505	Left Female	3"	76mm	182mm	150mm	70mm	510mm
RPL 506	Right Female	3"	76mm	182mm	150mm	70mm	510mm
RPL 505M	Left Male	3"	76mm	182mm	150mm	69mm	490mm
RPL 506M	Right Male	3"	76mm	182mm	150mm	69mm	490mm
RPL 507	Left Female	2¾"	70mm	168mm	135mm	70mm	500mm
RPL 508	Right Female	2¾"	70mm	168mm	135mm	70mm	500mm
RPL 511	Left Child	2½"	64mm	152mm	128mm	56mm	420mm
RPL 512	Right Child	2½"	64mm	152mm	128mm	56mm	420mm
RPL 513	Left Child	2¼"	57mm	132mm	126mm	48mm	340mm
RPL 514	Right Child	2¼"	57mm	132mm	126mm	48mm	340mm

Realistic Prosthetics Ltd

Unit 19 Platts Eyot, Lower Sunbury Road,
Hampton, Middlesex TW12 2HF UK

Phone +44 (0) 20 8783 0099

Fax +44 (0) 20 8941 4077

E-mail info@realisticprosthetics.com
realisticprosthetics.com

These cosmetic gloves are compatible with UK,
European and other manufacturers hand systems

CE marked product

For sales, please contact us directly,
or our distribution Agent in your Country.